

## DMEM/F12 with HEPES and L-glutamine

### High quality cell culture performance

Captivate Bio's DMEM/F12 Medium with HEPES is formulated for superior quality, reliability, and consistency for use in mammalian cell cultures. DMEM/F12 is a widely used serum-free, basal medium for supporting the growth of various stem cell cultures including human pluripotent, stem cells, induced pluripotent stem cells, human endothelial cells, and mouse embryonic stem cells.

#### This DMEM/F12 is modified as follows:

##### With

- HEPES
- L-glutamine
- Sodium Pyruvate

##### Without


- Phenol Red

#### Storage and handling

We recommend that this media be stored at 2-8°C and protected from light. This formulation uses a sodium bicarbonate buffer system and will require a 5-10% CO<sub>2</sub> environment to maintain physiological pH. Always use aseptic techniques when handling and supplementing cell culture media.

Product Information	CAT NO	QTY
DMEM/F12 Medium	DFHG01	500 mL
	DFHG02	1000 mL

\* Custom formulations and bulk packaging available.



High Quality Medium



Made in the USA



Reliable Formulation



Customizable

#### Quality Specifications

<b>Intended use</b>	Mammalian cell culture, hPSCs, fibroblasts, human, mouse, rat, non-human primate, other cells.		
<b>Glutamine</b>	L-glutamine	<b>HEPES Buffer</b>	HEPES
<b>Phenol Red Indicator</b>	No Phenol Red	<b>Shelf-Life</b>	12 months from date of manufacture
<b>pH</b>	7.0-7.4	<b>Storage</b>	2-8° C
<b>Osmolality</b>	280-305 mOsm/kg	<b>Sterility</b>	Pass
<b>Endotoxin</b>	< 1.0 EU/mL	<b>Mycoplasma</b>	Not Detected

For more information, visit [captivatebio.com](https://www.captivatebio.com), email [orders@captivatebio.com](mailto:orders@captivatebio.com), or contact us at (617) 607-4017.

## DMEM/F12 with HEPES, with L-glutamine, without Phenol Red

Catalog Number(s): DFHG01, DFHG02

Components	Molecular Weight	mg/L	mM
<b>Amino Acids</b>			
Glycine	75	18.75	0.25
L-Alanine	89	4.45	0.05
L-Arginine hydrochloride	211	147.5	0.69905216
L-Asparagine-H2O	150	7.5	0.05
L-Aspartic acid	133	6.65	0.05
L-Cystine hydrochloride-H2O	176	17.56	0.09977272
L-Cystine 2HCl	313	31.29	0.09996805
L-Glutamic Acid	147	7.35	0.05
L-Glutamine	146	365	2.5
L-Histidine hydrochloride-H2O	210	31.48	0.14990476
L-Isoleucine	131	54.47	0.41580153
L-Leucine	131	59.05	0.45076334
L-Lysine hydrochloride	183	91.25	0.4986339
L-Methionine	149	17.24	0.11570469
L-Phenylalanine	165	35.48	0.2150303
L-Proline	115	17.25	0.15
L-Serine	105	26.25	0.25
L-Threonine	119	53.45	0.44915968
L-Tryptophan	204	9.02	0.04421569
L-Tyrosine disodium salt dihydrate	261	55.79	0.21375479
L-Valine	117	52.85	0.4517094
<b>Vitamins</b>			
Biotin	244	0.0035	1.43E-05
Choline chloride	140	8.98	0.06414285
Calcium D-pantothenate	477	2.24	0.00469602
Folic Acid	441	2.65	0.00600907
Niacinamide	122	2.02	0.01655738
Pyridoxine hydrochloride	206	2	0.00970874
Riboflavin	376	0.219	5.82E-04
Thiamine hydrochloride	337	2.17	0.00643917
Vitamin B12	1355	0.68	5.02E-04
i-Inositol	180	12.6	0.07

Components	Molecular Weight	mg/L	mM
<b>Inorganic Salts</b>			
Calcium Chloride (CaCl <sub>2</sub> ) (anhyd.)	111	116.6	1.0504504
Cupric sulfate (CuSO <sub>4</sub> ·5H <sub>2</sub> O)	250	0.0013	5.20E-06
Ferric Nitrate (Fe(NO <sub>3</sub> ) <sub>3</sub> ·9H <sub>2</sub> O)	404	0.05	1.24E-04
Ferric sulfate (FeSO <sub>4</sub> ·7H <sub>2</sub> O)	278	0.417	0.0015
Magnesium Chloride (anhydrous)	95	28.64	0.30147368
Magnesium Sulfate (MgSO <sub>4</sub> ) (anhyd.)	120	48.84	0.407
Potassium Chloride (KCl)	75	311.8	4.1573334
Sodium Bicarbonate (NaHCO <sub>3</sub> )	84	1200	14.285714
Sodium Chloride (NaCl)	58	6995.5	120.61207
Sodium Phosphate dibasic (Na <sub>2</sub> HPO <sub>4</sub> ) anhydrous	142	71.02	0.50014085
Sodium Phosphate monobasic (NaH <sub>2</sub> PO <sub>4</sub> ·H <sub>2</sub> O)	138	62.5	0.45289856
Zinc sulfate (ZnSO <sub>4</sub> ·7H <sub>2</sub> O)	288	0.432	0.0015
<b>Other Components</b>			
D-Glucose (Dextrose)	180	3151	17.505556
HEPES	238	3574.5	15.018908
Hypoxanthine Na	159	2.39	0.01503145
Linoleic Acid	280	0.042	1.50E-04
Lipoic Acid	206	0.105	5.10E-04
Putrescine 2HCl	161	0.081	5.03E-04
Sodium Pyruvate	110	55	0.5
Thymidine	242	0.365	0.00150826