

Cytogenetics

High quality cell culture media and supplements for successful analysis

Captivate Bio offers a specialized portfolio of cell culture media, salts, buffers, FBS, and reagents for cytogenetic analysis. Products include BIO-AMF™, BIO-PB™, and BIO-HEMATO™ cell culture media which is optimized to accelerate cell growth with high mitotic index and superior morphology.

Each media is designed to produce clear results in the analysis of hematopoietic cells, peripheral blood, bone marrow, amniotic fluid, and chorionic villi (CVS) samples.



Complete, ready-to-use media



Save time with fast growth and rapid analysis



cGMP manufactured and performance tested



Quality control tested for sterility, pH, osmolality, and endotoxin

Product Information	CAT NO	QTY
BIO-AMF™-2 Complete Medium	01-194-1A	500 mL
	01-194-1B	100 mL
BIO-HEMATO™ Karyotyping Medium	01-200-1A	500 mL
	01-200-1B	100 mL
BIO-PB™ Karyotyping Medium with PHA-M	01-201-1A	500 mL
	01-201-1B	100 mL
Colcemid Solution	12-004-1D	10 mL
Potassium Chloride (KCl) Solution (0.075 molar)	12-005-1A	500 mL
Gentamicin Sulfate Solution	03-035-1B	100 mL
Cell Synchronization Kit	12-008-60	1 Kit
RPMI 1640, without glucose	01-101-1A	500 mL
RPMI 1640, with L-glutamine	RPMI01	500 mL
HBSS, with calcium and magnesium	HBS002	1 L
DPBS, without calcium and magnesium	DPBS02	1 L
FBS, US Origin	S1480	500 mL
FBS, USDA Origin	S1620	500 mL

Additional size and formats available on captivatebio.com.

For more information, visit captivatebio.com, email orders@captivatebio.com, or contact us at (617) 607-4017.