

# VTM Collection Kits

## Innovative technology for sample collection

Proper specimen collection and transport plays a critical role in laboratory diagnosis of infectious diseases associated with viruses. The Ruhof VTM Collection Kits offered by Captivate Bio allow for the collection, preservation and transport of viruses, such as coronavirus (COVID-19).

Captivate Bio is an authorized distributor of the Ruhof Corporation's viral transport medium and flocked swab collection kit.

Ordering Information	Cat. No.	Size
VTM Collection Kits	345TK2	25 Kits
	345TK2-4CS	100 Kits

Please inquire for bulk pricing and larger case packs.

For additional product or technical information, please visit [www.captivatebio.com](http://www.captivatebio.com), email [orders@captivatebio.com](mailto:orders@captivatebio.com), or call customer service at (617) 607-4017.



Active Transport Medium



Made in the USA



Ergonomic Swab Design



Suited for PCR and EIA Assays

Quality Specifications			
Intended use	Collection and transport of biological specimens		
Tube dimensions	15 x 95 mm	Fill volume	3 mL
Tube size	10 mL	Shelf-Life	18 months from date of manufacture
pH	7.1 to 7.5	Storage	15 - 30°C
Swabs (fiber)	Sterilized and free of DNase, RNase, human DNA		
Swabs (shaft)	Nasopharyngeal flocked swab breakpoint at 100 mm		
Formulation (VTM)	Hanks' Balanced Salt Solution Heat-Inactivated FBS Gentamicin Sulfate, Amphotericin B Phenol Red		