

HBSS, with Calcium and Magnesium

Catalog Number(s): HBS001, HBS002

Components	mg/L
Inorganic Salts	
Calcium Chloride	140
Potassium Chloride	400
Potassium Phosphate, monobasic	60
Magnesium Chloride	100
Magnesium Sulfate	100
Sodium Chloride	8000
Sodium Phosphate, dibasic	48
Sodium Bicarbonate	350

Components	mg/L
Other Components	
D-Glucose	1000

HBSS, with Calcium, Magnesium, and Phenol Red

Catalog Number(s): HBSPH1, HBSPH2

Components	mg/L
Inorganic Salts	
Calcium Chloride	140
Potassium Chloride	400
Potassium Phosphate, monobasic	60
Magnesium Chloride	100
Magnesium Sulfate	100
Sodium Chloride	8000
Sodium Phosphate, dibasic	48
Sodium Bicarbonate	350

Components	mg/L
Other Components	
D-Glucose	1000
Phenol Red	10